

OCTOPUS[®] TE

Tissue Stabilizer



Significant Step Closer to Minimally Invasive Endoscopic CABG

Octopus® TE Tissue Stabilizer is the latest innovation from Intuitive Surgical, the market leader in minimally invasive beating heart surgery products.

Whether you are performing coronary procedures with a robotic approach or using a small thoracotomy, the Octopus® TE stabilizer is designed to facilitate your procedure.

The Octopus® TE device is a reusable tissue stabilizer with multiple flexible pods that enable insertion into and removal from the thoracic cavity through a port, thus eliminating the need for an incision for insertion of the stabilizer.

Custom position

Stabilize



Collapsible pods for insertion/removal through a 12mm port

ORDERING INFORMATION

Product Code	Description	U/M
TS502	Octopus® TE Tissue Stabilizer	1 ea - Reusable
MR103	Universal Mounting Rail (Stabilizer)	1 ea - Reusable
DS501	Octopus® TE Irrigator	5/ctn
DS502	Octopus® TE Vacuum Lines	5/ctn
DS503	Octopus® TE Articulating Arm	5/ctn

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

For more information contact your local Medtronic Revascularization and Surgical Therapies Sales Representative or call Customer Service toll-free at 1-800-328-1357.

World Headquarters

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel: (763) 514-4000
Fax: (763) 514-4879
www.medtronic.com

Medtronic USA, Inc.
Toll-free: 1 (800) 328-2518
(24-hour technical support for
physicians and medical professionals)

Europe

Medtronic Europe Sàrl
Route du Molliau 31
CH-1131 Tolochenaz
Switzerland
Tel: (41 21) 802 7000
Fax: (41 21) 802 7900
www.medtronic.com

Canada

Medtronic of Canada Ltd.
6733 Kitimat Road
Mississauga, Ontario L5N 1W3
Canada
Tel: (905) 826-6020
Fax: (905) 826-6620
Toll-free: 1 (800) 268-5346

Asia Pacific

Medtronic International, Ltd.
16/F Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel: (852) 2891 4456
Fax: (852) 2891 6830
enquiryap@medtronic.com
www.medtronic.com

Latin America

Medtronic USA, Inc.
Doral Corporate Center II
3750 NW 87th Avenue Suite 700
Miami, FL 33178
USA
Tel: (305) 500-9328
Fax: (786) 709-4244
www.medtronic.com

UC200401703b EN
© Medtronic, Inc. 2005
All Rights Reserved
Printed in USA
August 2005

